

(425)355-2034 / 1(800)355-2036

www.nafzigmicroscopes.com

Meiji Microscopes

Discontinued

“FLUORESCENT” MODEL TM-420F

Discontinued

COMPOUND MICROSCOPE MODEL TM-208



- INCLINED 30 degrees monocular head with 360 degrees rotation and multicoated high transmission prism.
- Widefield WF 10x eyepiece with pointer 15.5mm eye-point-locked in.
- DIN Achromatic objectives with retractable mount on 20X and above.
- Built-in 0.65 NA condenser with disc diaphragm or 1.25 NA Abbe condenser with iris diaphragm on spiral mount.
- Large stand 119mm x 124mm with locked in stage clips.
- Rapid action coarse focus with rack and pinion and separate lever type fine focus.
- Auto focus stop and slip clutch to prevent rack and pinion damage.
- Coarse focus tension adjustment.
- Triple or quadruple objective changer.
- Rugged, well balanced stand 150mm(W), 210mm(D), 320mm(H), with non-slip rubber feet.

- Built-in Super lite illuminator with extra long life bulb. (up to 10,000 hrs.) and a “cool” minimum heat build-up, with approved 115V operation.
- Chemical resistant light grey leathertone finish.



- 10X wide field eyepiece locked-in with pointer.
- 40X objective spring retractable to prevent slide damage.
- Built-in NA 0.65 condenser with disc diaphragm.
- Built-in focus slip-clutch to prevent rack and pinion damage.
- Tension adjustment collar on course focus.
- Triple-objective changer.
- Bolted-on attachable illuminator with 20W bulb.
- Locked-in stage clips.
- Supplied with dust cover and manual.

1-4.....\$425.00

5 OR MORE.....\$395.00

(LIST PRICE \$495.00)

1-4.....\$415.00

5 OR MORE.....\$385.00

(LIST PRICE \$480.00)

STEREO MICROSCOPE MODEL SKT SERIES



- Crisp, erect images over a wide field of view.
- Binocular body is inclined at 45 degrees for great ease of observation.
- Binocular body can be raised or lowered on pillar stand.
- Interpupillary distance is adjustable between 51mm and 75mm.
- Dioptic adjustment is possible with one eyepiece tube.
- Focusing is by rack and pinion with safety stopper at the lower end of traverse.
- Equipped with pairs of objectives TURRET MOUNT for rapid magnification range.
- Magnification range from 5X to 80X with optional eyepieces.
- Stages are equipped with stage clips.
- Illumination employs 12V, 10W bulb.

- Incident and transmitted lights are designed to operate separately.
- Supplied in styrofoam.

SKT-1BT--1X,2X Objectives

SKT-2BT--1X,3X Objectives

SKT-3BT--2X,4X Objectives

1-4.....\$355.00

5 OR MORE.....\$335.00

(LIST PRICE \$450.00)

STEREO MICROSCOPE MODEL EMT-3PB



- Crisp, erect image, with flat, widefield of view and long working distance.
- Paired widefield objectives in a turret mount for rapid selection of magnification (EMT).
- Binocular tubes with an inclination of 45 degrees for ease of observation.
- Interpupillary distance adjustment from 54mm to 75mm.
- Dioptic adjustment is possible on the left eyepiece tube.
- Binocular body can be rotated a full 360 degrees, and locked in any position.

1-4.....\$725.00

5 OR MORE.....\$695.00

(LIST PRICE \$950.00)