

(425)355-2034 / 1(800)355-2036

www.nafzigmicroscopes.com

National Magnavision

MICROSCOPES WITH BUILT-IN DIGITAL CAMERA



MODEL DC3-163



FIVE EXCITING NEW MICROSCOPES WITH BUILT-IN DIGITAL CAMERAS! EACH MODEL IS A COMPLETE PACKAGE, READY TO GO AND EASY TO USE

More than just great microscopes.....simply plug them into your computer and you are ready to capture images of your specimen, label them, measure them, store them or e-mail them. You don't need to buy a separate microscope, camera, cables, and software. Everything needed to get started is contained within each package.

All models include the same easy-to-use software for Windows 98 and feature the same computer applications. You can also record audio comments along with your images. Two models also permit simultaneous hook up to video monitor and to a large screen projector. See individual model descriptions for details. USB, S and RCA cables, as required for each individual model, are included along with complete instructions for easy "plug and play" use.

DC SERIES MICROSCOPES WITH BUILT-IN DIGITAL CAMERAS

Prices include microscope with built in digital camera, software CD for Windows 98 or above and all cables required for operation.

Model No.	Description	Bulb	No. List Price	Your Price
DC4-155	Compound microscope, 4x, 10x, 40xR objectives, separate coarse & fine focus, 12v 20 watt halogen illumination, 0.65 in-stage condenser, digital signal only.	800-160	\$ 749.00	\$ 550.00
DC4-156-S	Compound microscope, 4x, 10x, 40xR, 100xR objectives, separate coarse & fine focus, 12v 20 watt halogen illumination, 1.25 N.A. Abbe condenser, digital signal only.	800-160	\$ 949.00	\$ 750.00
DC3-163	Compound microscope, 4x, 10x, 40xR, 100xR objectives, coaxial coarse & fine focus, 12v 20 watt halogen illumination, 1.25 N.A. Abbe condenser, digital & analog signals.	800-160	\$1499.00	\$1190.00
DC4-456H	Stereo microscope, 10x & 30x magnification, 12v 10 watt tungsten illumination, digital signal only.	800-400	\$ 749.00	\$ 550.00
DC3-420TH	Stereo zoom microscope, 10x - 40x magnification, top illumination-12v15 watt halogen, bottom illumination-12v 10 watt halogen, digital & analog signals.	800-421(T) 800-170(B)	\$1475.00	\$1180.00