

(425)355-2034 / 1(800)355-2036

www.nafzigmicroscopes.com

National Video Cameras

MOTICAM 352 DIGITAL CAMERA with Motic Images 2000 Software

- Digital camera, with USB cable attached, 1/4 in. chip
- Digital signal, 3 Mega Pixels
- C-ring camera adapter
- CCD 8mm lens
- 4 eyepiece adapters, fit O.D. eyepieces, 28mm, 30mm, 34mm, 35mm
 - 2 calibration slides (compound and stereo)
- Motic images 2000 software
- Motic images 2000 J 1.0 software for Mac OS 9.2



OUR PRICE.....\$170.00

(List Price \$199.00)

DAC-3 DIGITAL/ANALOG CAMERA with Motic Images 2000 Software

- Digital/analog camera, 1/3 in.
- Analog signal 480 lines of TV resolution
- Digital signal, 4 Mega Pixels
- C-ring adapter
- 12mm CCD lens
- 4 eyepiece adapters to fit eyepiece O.D.'s of 28mm, 30mm, 34mm, 35mm
- USB, RCA, and S-video cables
 - 110v-240v power converter and cord
- 2 calibration slides (compound and stereo)
- Motic images 2000 Plus software
- Motic images 2000 J1.0 software for Mac OS 9.2



OUR PRICE.....\$599.00

(List Price \$699.00)

