

(425)355-2034 / 1(800)355-2036

www.nafzigmicroscopes.com

# National Video Cameras

## D-MOTICAM 483



- Digital / Analog Camera 1/3 inch
- Analog Signal (S-Video, RCA) - 640 x 480 resolution
- Digital Signal (USB) - 640 X 480 resolution
- Minimum illumination: 3 Lux (F.2)
- Data and power sockets are integrated in the unit base
- 12 mm CCD Lens
  - 28 mm, 30mm, 34mm, 35mm, eyepiece adapters included
- USB, RCA and S-Video cables included
- Switching 110v - 240v power converter and cord
- Hardware white balance adjustments
- C-Ring Adapter
- 2 calibration slides
- Motic Software CD
- Carrying / Storage Case

**OUR PRICE.....\$599.00**

(List Price \$749.00)



**(Laptop Not Included)**