

(425)355-2034 / 1(800)355-2036

www.nafzigmicroscopes.com

Swift Microscopes

“Fluorescent” Compound Microscope MODEL M2251F



- Widefield eyepiece, W10X, with built-in pointer, locked-on.
- Inclined eye tube of constant height.
- Objective, 4X, 10X, 40XRD; Triple nosepiece.
- Monocular body is 360 degrees rotatable.
- Large stage, 115mm x 120mm.
- 5w Fluorescent Bulb
- Locked-on stage clips.
- Disc diaphragm with five apertures.
- Separate coarse and fine focus controls.
- Tension control clutch system.
- Approved three-wire grounded electrical system.
- Limited lifetime warranty.

1-4.....\$220.00
5 OR MORE.....\$200.00
(LIST PRICE \$249.00)

NEW

“Cordless” Compound Microscope MODEL M2251C



- Up to 50 hours of continuous use on a single, 8-hour charge
- Battery life of 500 recharge cycles (approximately 2.4 years)
- Energy-efficient technology produces up to 100,000 hours of light equal to that of a 20 watt bulb
- Dual-purpose AC adapter/recharger (included) facilitates uninterrupted function
- M2251C can be used while battery is recharging
- Portability of M2251C makes it suitable for both classroom use and off-site or outdoor study

- Locked-on, W10X eyepiece and stage clips discourage theft
- Separate coarse and fine focus system is easy to use and simple to repair
- Slip-clutch focus system protects gears from over-focusing
- Retractable 40XRD objective protects delicate lenses from impact damage
- The M2251C is covered by Swift's unequalled Lifetime Warranty (excluding battery and bulb)

1-4.....\$240.00
5 OR MORE.....\$220.00
(LIST PRICE \$269.00)

“Fluorescent” Compound MODEL M3501DF



- One-piece eyetube and condenser housing.
- Plastic shield for iris diaphragm protection.
- Locked-on eyetube and stage clips.
- Cord-hanger allows for safe and hassle-free cord storage.
- Reversed nosepiece allows for easy slide manipulation.
- MA889A (high-drive) and MA916 (low-drive) mechanical stages (optional).
- Phase Contrast “Quodlibet” condenser

- and objectives (optional).
- Built to withstand years of wear and tear!
- Eternal parts support.
- Limited Lifetime Warranty on all parts and workmanship.
- UL/CUL-approved electrical system.
- Electrical cord is grommeted to the base to avoid cord detachment.
- UltraLite the original and brightest “cool light” illuminator.
- Pointmaster eyepiece for measuring specimens.
- Unique “C-wrench” keeps control of focusing tension in teacher’s hands.
- Savings on servicing and parts replacement common to lesser-quality brands.
- Ultra-Lite illumination reduces your electrical and bulb replacement costs.

1-4.....\$295.00
5 OR MORE.....\$280.00
(LIST PRICE \$349.00)

SWIFT FIVE

Advanced Binocular Microscopes



- Siedentopf Head provides larger interpupillary range and is ergonomically positioned for operator comfort
- Focus-tension control knob offers “feathertouch” sensitivity for extremely precise focusing
- All-metal internal parts and high-grade synthetic lubricants provide a smooth, precise operation.
- Heads are inclined at 45 degrees and rotatable 360 degrees for ergonomic viewing
- Inward-facing (reversed) nosepiece allows user clear access to stage for

- slide manipulation.
- 40XRD & 100XRD objectives are spring-loaded and sealed to protect optics from contact damage.
- Large, research-level, mechanical stage accommodates two (2) slides at a time for comparison viewing
- Halogen illumination with variable intensity control provides optimal light

List Price Our Price

N5S..(Semi-Plan Obj).....\$1,170.....\$960.00
M5P..(Plan-Achro Obj).....\$1,295.....\$1040.00